

FAM110D RABBIT PAB

Cat.#: S215923

Product Name: Anti-FAM110D Rabbit Polyclonal Antibody

Synonyms: GRRP1

UNIPROT ID: Q8TAY7 (Gene Accession - NP_079145)

Background: The function of this protein remains unknown.

Immunogen: Synthetic peptide of human FAM110D

Applications: ELISA, IHC

Recommended Dilutions: IHC: 30-150; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human tonsil tissue using 215923(FAM110D Antibody) at a dilution of 1/60(Cytoplasm and Cell membrane).

In comparison with the IHC on the left, the same paraffin-embedded Human tonsil tissue is first treated with the synthetic peptide and then with 215923(Anti-FAM110D Antibody) at dilution 1/60.