

FAM171A2 RABBIT PAB

Cat.#: S222167

Product Name: Anti-FAM171A2 Rabbit Polyclonal Antibody

Synonyms:

UNIPROT ID: A8MVW0 (Gene Accession - NP_940877)

Background: The function of this protein remains unknown.

Immunogen: Synthetic peptide of human FAM171A2

Applications: ELISA, IHC

Recommended Dilutions: IHC: 20-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human tonsil tissue using 222167(FAM171A2 Antibody) at a dilution of 1/40(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human tonsil tissue is first treated with the synthetic peptide and then with 222167(Anti-FAM171A2 Antibody) at dilution 1/40.