

FAM98A RABBIT PAB

Cat.#: S215780

Product Name: Anti-FAM98A Rabbit Polyclonal Antibody

Synonyms:

UNIPROT ID: Q8NCA5 (Gene Accession - NP_056290)

Background: FAM98A belongs to the FAM98 family. The function of FAM98A remains unknown.

Immunogen: Synthetic peptide of human FAM98A

Applications: ELISA, WB, IHC

Recommended Dilutions: IHC: 30-150;WB: 200-1000;ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human thyroid cancer tissue using 215780 (FAM98A Antibody) at a dilution of 1/20 (Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human thyroid cancer tissue is first treated with the synthetic peptide and then with 215780 (Anti-FAM98A Antibody) at dilution 1/20.



Gel: 8% SDS-PAGE, Lysate: 40 μ g;
Lane 1-3: Human breast cancer tissue, Human liver cancer tissue and Human lung adenocarcinoma tissue lysates;
Primary antibody: 215780 (FAM98A Antibody) at dilution 1/300;
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution;
Exposure time: 10 minutes