

FAS RABBIT PAB

Cat.#: N225763

Product Name: Anti-Fas Rabbit pAb

Synonyms: FAS; ALPS1A; APO-1; APT1; CD95; FAS1; FASTM; TNFRSF6

UNIPROT ID: P25445

Background: Receptor for TNFSF6/FASLG. The adapter molecule FADD recruits caspase-8 to the activated receptor. The resulting death-inducing signaling complex (DISC) performs caspase-8 proteolytic activation which initiates the subsequent cascade of caspases (aspartate-specific cysteine proteases) mediating apoptosis.

Immunogen: A synthesized peptide derived from human Fas

Applications: WB,IHC-P,ICC/IF,FC

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 38 kDa; Observed MW: 48 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Fas in Ramos lysates using Fas antibody.