

FEN1 (7H8) MOUSE MAB

Cat.#: N261126

Product Name: Anti-FEN1 (7H8) Mouse Monoclonal Antibody

Synonyms: FEN1; RAD2; Flap endonuclease 1; FEN-1; DNase IV; Flap structure-specific endonuclease 1; Maturation factor 1; MF1; hFEN-1

UNIPROT ID: P39748

Background: Acts as a genome stabilization factor that prevents flaps from equilibrating into structures that lead to duplications and deletions. Also possesses 5'-3' exonuclease activity on nicked or gapped double-stranded DNA, and exhibits RNase H activity. Also involved in replication and repair of rDNA and in repairing mitochondrial DNA.

Immunogen: Purified recombinant human FEN-1 protein fragments expressed in E.coli.

Applications: WB, ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IF: 1/50-1/200

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 7H8-F4-C11

MW: Calculated MW: 43 kDa; Observed MW: 43 kDa

Isotype: IgG1

Purification: Affinity Purified

Species Reactivity: Human, Mouse, Rat, Monkey, Hamster


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of FEN1 (7H8) in HeLa using FEN1 antibody.



Western blot analysis of FEN1 in HeLa, Jurkat, 3T3, COS7, PC-12, C6, Raw264.7, Raji and CHO-K1 lysates using FEN1 antibody.