

FETUIN A RABBIT MAB

Cat.#: N262216

Product Name: Anti-Fetuin A Rabbit Monoclonal Antibody

Synonyms: 16G2; A2HS; AHS; Ahsg; FETUA; FETUB; FetuinA; FetuinB; Gugu; HSGA; IRL685; PRO2743

UNIPROT ID: P02765

Background: Promotes endocytosis, possesses opsonic properties and influences the mineral phase of bone. Shows affinity for calcium and barium ions.

Immunogen: A synthetic peptide of human Fetuin A

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R03-5A4

MW: Calculated MW: 39 kDa; Observed MW: 45-50 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Fetuin A in HepG2, 3T3, U937 lysates using Fetuin A antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
