

FGFR4 MOUSE MAB

Cat.#: N261191

Product Name: Anti-FGFR4 Mouse Monoclonal Antibody

Synonyms: FGFR4; JTK2; TKF; Fibroblast growth factor receptor 4; FGFR-4; CD334

UNIPROT ID: P22455

Background: Receptor for acidic fibroblast growth factor. Does not bind to basic fibroblast growth factor. Binds FGF19.

Immunogen: Purified recombinant extracellular fragment of human FGFR4 fused with hlgGfc tag expressed in HEK293 cell line.

Applications: ICC/IF

Recommended Dilutions: IF: 1/50-1/200

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 7H1

MW: -

Isotype: IgG1

Purification: Ascitic Fluid

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunofluorescence analysis of FGFR4 in HEK293 cells transfected with FGFR4hlgGfc using FGFR4 antibody (green) Blue: DRAQ5 fluorescent DNA dye.



Western blot analysis of FGFR4 in Human FGFR4 (aa22369) lysates using FGFR4 antibody.