

FKBP51 RABBIT PAB

Cat.#: N225978

Product Name: Anti-FKBP51 Rabbit pAb

Synonyms: P54; AIG6; FKBP51; FKBP54; PPIase; Ptg-10

UNIPROT ID: Q13451

Background: Interacts with functionally mature heterooligomeric progesterone receptor complexes along with HSP90 and TEBP.

Immunogen: A synthesized peptide derived from human FKBP51

Applications: WB,IHC-P,ICC/IF,IP,FC

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 51 kDa; Observed MW: 51 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FKBP51 in Caco2 lysates using FKBP51 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
