

FLAG TAG MAB**Flag Tag 小鼠单克隆抗体**

货号: 28002

规格: 100 μ L

Symbol: FLAG(M2)

描述: Flag Tag 小鼠单克隆抗体

背景: The Flag antibody recognizes any protein containing Flag Tag.

免疫原: DYKDDDDKC

经过测试的应用: ELISA, WB, IF, IHC

推荐稀释度:

ELISA 1:2000-1:5000

WB 1:1000-1:5000

IF 1:100-1:200

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG


纯化:: Purified from ascites

Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Recognizes protein contain Flag.

储存条件:: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:

Western blot: Recombinant Flag tag fusion protein was diluted with Flag mouse monoclonal antibody (1:1000), and the protein loading amount was 50ng, and Western blot analysis was performed.

