

FLT3 LIGAND RABBIT PAB

Cat.#: N225190

Product Name: Anti-Flt3 Ligand Rabbit pAb

Synonyms: FLT3LG; Fms-related tyrosine kinase 3 ligand; Flt3 ligand; Flt3L; SL cytokine

UNIPROT ID: P49771

Background: Stimulates the proliferation of early hematopoietic cells by activating FLT3.

Immunogen: The antiserum was produced against synthesized peptide derived from the C-terminal region of human FLT3LG. AA range:171-220

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 26 kDa; Observed MW: 26 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Flt3 Ligand in HeLa lysates using Flt3 Ligand antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
