

FOXP1 RABBIT MAB

Cat.#: N262236

Product Name: Anti-FOXP1 Rabbit Monoclonal Antibody

Synonyms: FOXP1; HSPC215; Forkhead box protein P1

UNIPROT ID: Q9H334

Background: Transcriptional repressor. It plays an important role in the specification and differentiation of lung epithelium. Can act with CTBP1 to synergistically repress transcription but CTBP1 is not essential. Essential transcriptional regulator of B cell development.

Immunogen: A synthetic peptide of human FOXP1

Applications: WB,IHC-F,IHC-P,ICC/IF,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R01-6I3

MW: Calculated MW: 75 kDa; Observed MW: 82-90 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of FOXP1 in Raw264.7, rat lung lysates using FOXP1 antibody.