

GAS7 RABBIT PAB

Cat.#: S210332

Product Name: Anti-GAS7 Rabbit Polyclonal Antibody

Synonyms: MLL/GAS7

UNIPROT ID: O60861 (Gene Accession - BC001152)

Background: Growth arrest-specific 7 is expressed primarily in terminally differentiated brain cells and predominantly in mature cerebellar Purkinje neurons. GAS7 plays a putative role in neuronal development. Several transcript variants encoding proteins which vary in the N-terminus have been described.

Immunogen: Fusion protein of human GAS7

Applications: ELISA, WB, IHC

Recommended Dilutions: IHC: 50-200;WB: 500-2000;ELISA: 2000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling, Cancer


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human thyroid cancer tissue using 210332(GAS7 Antibody) at a dilution of 1/50(Cytoplasm).



Immunohistochemistry analysis of paraffin-embedded Human brain tissue using 210332(Anti-GAS7 Antibody) at a dilution of 1/50.



Gel: 8%SDS-PAGE, Lysate: 40 μ g;
Lane 1-2: Human fetal brain tissue, mouse brain tissue;
Primary antibody: 210332(GAS7 Antibody) at dilution 1/625;
Secondary antibody: Goat anti rabbit IgG at 1/8000 dilution;
Exposure time: 2 minutes