

**GAZ PAB****Gα<sub>z</sub> Polyclonal Antibody**

货号: 21016

基因符号: GNAZ

描述: Gα<sub>z</sub> Rabbit Polyclonal Antibody

背景: Heterotrimeric G proteins are critical cellular signal transducers. Gα<sub>z</sub> belongs to the Gi sub-family of G proteins that could mediate the inhibition of adenylyl cyclases. Other biochemical and physiological functions of Gα<sub>z</sub> proteins are being explored.

免疫原: Recombinant full length protein of Gα<sub>z</sub>.

经过测试的应用: ELISA, WB, IHC

推荐稀释度: ELISA 1:100-1:5000; WB 1:100-1:2000

浓度: 1 mg/ml

种属反应性: Rabbit

形式: Liquid

克隆性: Polyclonal

亚型: IgG


纯化: Purified from serum

**Preservative:** No

成分: PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Gα<sub>z</sub> antibody recognizes Gα<sub>z</sub> of vertebrates.

储存条件: Store at -20°C. Avoid repeated freezing and thawing

**Western blot:****kDa****WB: anti-G $\alpha_z$  pAb**

Western blot analysis of G $_z$  proteins. HEK-293 whole-cell lysate was used for western blots.