

GLYAT RABBIT PAB

Cat.#: S218814

Product Name: Anti-GLYAT Rabbit Polyclonal Antibody

Synonyms: CAT; GAT; ACGNAT

UNIPROT ID: Q6IB77 (Gene Accession - BC009785)

Background: The glycine-N-acyltransferase protein conjugates glycine with acyl-CoA substrates in the mitochondria. The protein is thought to be important in the detoxification of endogenous and xenobiotic acyl-CoA's. Two transcript variants encoding different isoforms have been found for this gene.

Immunogen: Fusion protein of human GLYAT

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Metabolism, Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human liver cancer tissue using 218814 (GLYAT Antibody) at a dilution of 1/30 (Cytoplasm).