

GSDME RABBIT MAB

Cat.#: N262128

Product Name: Anti-GSDME Rabbit Monoclonal Antibody

Synonyms: DFNA5; ICERE-1; GSDME

UNIPROT ID: O60443

Background: Plays a role in the TP53-regulated cellular response to DNA damage probably by cooperating with TP53 (PubMed:16897187, PubMed:18223688).

Immunogen: Recombinant protein of human DFNA5/GSDME

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-511

MW: Calculated MW: 55 kDa; Observed MW: 55 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of DFNA5/GSDME in U251MG, U2OS, HepG2, HeLa lysates using GSDME antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
