

GST3 RABBIT PAB

Cat.#: N225988

Product Name: Anti-GST3 Rabbit pAb

Synonyms: GSTP1; FAES3; GST3; Glutathione S-transferase P; GST class-pi; GSTP1-1

UNIPROT ID: P09211

Background: Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Immunogen: A synthesized peptide derived from human GST3 / GST pi

Applications: WB,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 23 kDa; Observed MW: 23 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of GST3 / GST pi in PC-3 lysates using GST3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
