

## GSTCD RABBIT PAB

**Cat.#:** S219163

**Product Name:** Anti-GSTCD Rabbit Polyclonal Antibody

**Synonyms:**

**UNIPROT ID:** Q8NEC7 (Gene Accession - BC032942 )

**Background:** The function of this protein remains unknown.

**Immunogen:** Fusion protein of human GSTCD

**Applications:** ELISA, WB, IHC

**Recommended Dilutions:** IHC: 70-350;WB: 500-2000;ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Cell Biology


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 219163(GSTCD Antibody) at a dilution of 1/70(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the fusion protein and then with 219163(Anti-GSTCD Antibody) at dilution 1/70.



Gel: 8%SDS-PAGE, Lysate: 40  $\mu$ g;  
Lane: K562 cell lysate;  
Primary antibody: 219163(GSTCD Antibody) at dilution 1/800;  
Secondary antibody: HRP-conjugated Goat anti rabbit IgG at 1/5000 dilution;  
Exposure time: 5 minutes