

GTPASE NRAS RABBIT MAB

Cat.#: N262816

Product Name: Anti-GTPase NRAS Rabbit Monoclonal Antibody

Synonyms: NS6; CMNS; NCMS; ALPS4; N-ras; NRAS1

UNIPROT ID: P01111

Background: Ras proteins bind GDP/GTP and possess intrinsic GTPase activity. Miscellaneous mutations which change AA 12, 13 or 61 activate the potential of Ras to transform cultured cells and are implicated in a variety of human tumors.

Immunogen: A synthetic peptide of human Ras

Applications: WB, IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-8E1

MW: Calculated MW: 21 kDa; Observed MW: 21 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human, Mouse, Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Ras in HeLa, 3T3, C6, rat Brain, MCF-7 lysates using GTPase NRAS antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
