

## GTPBP4 RABBIT MAB

**Cat.#:** N263353

**Product Name:** Anti-GTPBP4 Rabbit Monoclonal Antibody

**Synonyms:** CRFG; GTP binding protein 4; Gtpbp4; NGB; NOG1

**UNIPROT ID:** Q9BZE4

**Background:** Involved in the biogenesis of the 60S ribosomal subunit.

**Immunogen:** A synthetic peptide of human GTPBP4/NOG1

**Applications:** WB,IHC-P,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R07-8G0

**MW:** Calculated MW: 74 kDa; Observed MW: 80 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human tonsil using GTPBP4/NOG1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of GTPBP4/NOG1 in A549, HL-60 lysates using GTPBP4 antibody.