

HISTONE H1.2 RABBIT MAB

Cat.#: N262327

Product Name: Anti-Histone H1.2 Rabbit Monoclonal Antibody

Synonyms: H1.a; H1F2; H1s-1; HIST1H1C; Histone H1.2; Histone H1c; Histone H1d; Histone H1s-1

UNIPROT ID: P16403

Background: Histones H1 are necessary for the condensation of nucleosome chains into higher order structures.

Immunogen: A synthetic peptide of human Histone H1.2

Applications: WB,IHC-F,IHC-P,ICC/IF

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-6F1

MW: Calculated MW: 21 kDa; Observed MW: 30 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Histone H1.2 in HeLa, C6, MCF-7 lysates using Histone H1.2 antibody.



Western blot analysis of Histone H1.2 in mouse lung lysates using Histone H1.2 antibody.