

HISTONE H1.8 RABBIT PAB

Cat.#: N225481

Product Name: Anti-Histone H1.8 Rabbit pAb

Synonyms: H1.8; H1oo; osH1

UNIPROT ID: Q8IZA3

Background: May play a key role in the control of gene expression during oogenesis and early embryogenesis, presumably through the perturbation of chromatin structure. Essential for meiotic maturation of germinal vesicle-stage oocytes. The somatic type linker histone H1c is rapidly replaced by H1oo in a donor nucleus transplanted into an oocyte. The greater mobility of H1oo as compared to H1c may contribute to this rapid replacement and increased instability of the embryonic chromatin structure. The rapid replacement of H1c with H1oo may play an important role in nuclear remodeling .

Immunogen: The antiserum was produced against synthesized peptide derived from human Histone H1. AA range:131-180

Applications: IHC-P,IHC-F,ICC/IF,ELISA

Recommended Dilutions: IHC: 1/50-1/100 IF: 1/50-1/200 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: -

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunofluorescence analysis of Histone H1.8 in HeLa cells using Histone H1.8 antibody. Sample with blocking peptide on the right.