

HJURP RABBIT MAB

Cat.#: N263366

Product Name: Anti-HJURP Rabbit Monoclonal Antibody

Synonyms: FAKTS; URLC9; hFLEG1

UNIPROT ID: Q8NCD3

Background: Centromeric protein that plays a central role in the incorporation and maintenance of histone H3-like variant CENPA at centromeres.

Immunogen: Recombinant protein of human HJURP

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-4B7

MW: Calculated MW: 84 kDa; Observed MW: 84 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using HJURP antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of HJURP in HeLa, CHO-K1 lysates using HJURP antibody.