

HOXC9 RABBIT PAB

Cat.#: S221429

Product Name: Anti-HOXC9 Rabbit Polyclonal Antibody

Synonyms: HOX3; HOX3B

UNIPROT ID: P31274 (Gene Accession - NP_008828)

Background: This gene belongs to the homeobox family of genes. The homeobox genes encode a highly conserved family of transcription factors that play an important role in morphogenesis in all multicellular organisms. Mammals possess four similar homeobox gene clusters, HOXA, HOXB, HOXC and HOXD, which are located on different chromosomes and consist of 9 to 11 genes arranged in tandem. This gene is one of several homeobox HOXC genes located in a cluster on chromosome 12.

Immunogen: Synthetic peptide of human HOXC9

Applications: ELISA, WB, IHC

Recommended Dilutions: IHC: 50-200;WB: 200-1000;ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human breast cancer tissue using 221429(HOXC9 Antibody) at a dilution of 1/90(Nucleus).



Gel: 12%SDS-PAGE, Lysate: 40 µg;
Lane 1-4: C2C12, NIH/3T3, TM4 cell, Mouse skeletal muscle tissue lysates;
Primary antibody: 221429(HOXC9 Antibody) at dilution 1/350;
Secondary antibody: HRP-conjugated Goat anti rabbit IgG at 1/5000 dilution;
Exposure time: 30 seconds



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
