

## HSD17B1 RABBIT MAB

**Cat.#:** N262364

**Product Name:** Anti-HSD17B1 Rabbit Monoclonal Antibody

**Synonyms:** hydroxysteroid (17-beta) dehydrogenase 1; HSD17; EDHB17; EDH17B2; SDR28C1

**UNIPROT ID:** P14061

**Background:** Favors the reduction of estrogens and androgens. Also has 20-alpha-HSD activity. Uses preferentially NADH.

**Immunogen:** A synthetic peptide of human HSD17B1

**Applications:** WB,IHC-F,IHC-P,ICC/IF,IP

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20

**Host Species:** Rabbit

**Clonality:** Rabbit Monoclonal

**Clone ID:** R09-5K9

**MW:** Calculated MW: 35 kDa; Observed MW: 35 kDa

**Isotype:** IgG

**Purification:** Affinity Purified

**Species Reactivity:** Human


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of HSD17B1 in Human placenta lysates using HSD17B1 antibody.



# Product Description

Pioneering GTPase and Oncogene Product Development since 2010

---