



HSP60 RABBIT MAB**Cat.#:** N262369**Product Name:** Anti-Hsp60 Rabbit Monoclonal Antibody**Synonyms:** HSPD1; HSP60; 60 kDa heat shock protein; mitochondrial; 60 kDa chaperonin; Chaperonin 60; CPN60; Heat shock protein 60; HSP-60; Hsp60; HuCHA60; Mitochondrial matrix protein P1; P60 lymphocyte protein**UNIPROT ID:** P10809**Background:** HSP60 Implicated in mitochondrial protein import and macromolecular assembly. May facilitate the correct folding of imported proteins. May also prevent misfolding and promote the refolding and proper assembly of unfolded polypeptides generated under stress conditions in the mitochondrial matrix. Belongs to the chaperonin (HSP60) family. Interacts with HBV protein X and HTLV-1 protein p40tax.**Immunogen:** A synthetic peptide of human Hsp60**Applications:** WB,IHC-P**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100**Host Species:** Rabbit**Clonality:** Rabbit Monoclonal**Clone ID:** R02-5B7**MW:** Calculated MW: 61 kDa; Observed MW: 61 kDa**Isotype:** IgG**Purification:** Affinity Purified**Species Reactivity:** Human,Mouse,Rat**Conjugation:** Unconjugated**Modification:** Unmodified**Constituents:** PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide**Research Areas:** Tags & Cell Markers**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human lung cancer using Hsp60 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Hsp60 in 3T3, C6, mouse heart, rat heart lysates using Hsp60 antibody.



Western blot analysis of Hsp60 in Hela lysates using Hsp60 antibody.