

HUMAN SERUM ALBUMIN (5G11) MOUSE MAB

Cat.#: N261405

Product Name: Anti-Human Serum Albumin (5G11) Mouse Monoclonal Antibody

Synonyms: ALB; ALBU_HUMAN; Albumin (32 AA); Albumin (AA 34); Albumin; Analbuminemia; Bisalbuminemia; Cell growth inhibiting protein 42; DKFZp779N1935; Dysalbuminemic hyperthyroxinemia; Growth inhibiting protein 20; HSA; Hyperthyroxinemia dysalbuminemic; Serum albumin.

UNIPROT ID: P02768

Background: Albumin is a soluble, monomeric protein which comprises about one-half of the blood serum protein. Albumin functions primarily as a carrier protein for steroids, fatty acids, and thyroid hormones and plays a role in stabilizing extracellular fluid volume.

Immunogen: Native Human Serum Albumin..

Applications: LFIA,ELISA

Recommended Dilutions: ELISA: 1/10000

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 5G11-F5

MW: -

Isotype: IgG1

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Lateral Flow Immunoassay(LFIA) assay for antiHuman Serum Albumin mouse mAb (competitive assay).