

## IKAROS RABBIT PAB

**Cat.#:** N225942

**Product Name:** Anti-Ikaros Rabbit pAb

**Synonyms:** CLL associated antigen KW 6; DNA-binding protein Ikaros; hlk 1; hlk-1; Hs.54452; IK1; Ikaros (zinc finger protein); IKAROS; IKAROS family zinc finger 1 (Ikaros); Ikaros family zinc finger protein 1; Ikzf1; IKZF1\_HUMAN; LYF1; Lymphoid transcription factor LyF-1; PRO0758; Zinc finger protein subfamily 1A 1 (Ikaros); zinc finger protein subfamily 1A 1; Zinc finger protein; subfamily 1A; member 1; ZNFN1A1.

**UNIPROT ID:** Q13422

**Background:** Transcription regulator of hematopoietic cell differentiation (PubMed:17934067). Binds gamma-satellite DNA (PubMed:17135265, PubMed:19141594). Plays a role in the development of lymphocytes, B- and T-cells.

**Immunogen:** A synthesized peptide derived from human Ikaros

**Applications:** WB,IHC-P,FC

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 58 kDa; Observed MW: 50-70 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human,Mouse,Rat


**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Epigenetics and Nuclear Signaling

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of Ikaros in MOLT4 lysates using Ikaros antibody.