

IKZF4 RABBIT PAB

Cat.#: S221160

Product Name: Anti-IKZF4 Rabbit Polyclonal Antibody

Synonyms: EOS; ZNFNIA4

UNIPROT ID: Q9H2S9 (Gene Accession - NP_071910)

Background: Members of the Ikaros (ZNFNIA1; MIM 603023) family of transcription factors, which includes Eos, are expressed in lymphocytes and are implicated in the control of lymphoid development.

Immunogen: Synthetic peptide of human IKZF4

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 2000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 221160(IKZF4 Antibody) at a dilution of 1/25(Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the synthetic peptide and then with 221160(Anti-IKZF4 Antibody) at dilution 1/25.