

IL-17B RABBIT PAB

Cat.#: N225221

Product Name: Anti-IL-17B Rabbit pAb

Synonyms: IL-17B; interleukin 17B; NIRF; IL-20; ZCYTO7

UNIPROT ID: Q9UHF5

Background: Stimulates the release of tumor necrosis factor alpha and IL-1-beta from the monocytic cell line THP-1.

Immunogen: Synthetic peptide from human protein at AA range: 1-50

Applications: IHC-P,ELISA

Recommended Dilutions: IHC: 1/50-1/100 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: -

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human small intestine using IL-17B antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
