

IL-22 RABBIT PAB

Cat.#: N225206

Product Name: Anti-IL-22 Rabbit pAb

Synonyms: IL-22; IL22; ILTIF; ZCYTO18; Interleukin-22; IL-22; Cytokine Zcyto18; IL-10-related T-cell-derived-inducible factor; IL-TIF

UNIPROT ID: Q9GZX6

Background: Cytokine that contributes to the inflammatory response in vivo.

Immunogen: The antiserum was produced against synthesized peptide derived from the C-terminal region of human IL22. AA range:121-170

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 20 kDa; Observed MW: 20 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Immunology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of IL-22 in mouse KIDNEY lysates using IL-22 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
