

## INPP4B RABBIT PAB

**Cat.#:** N225707

**Product Name:** Anti-INPP4B Rabbit pAb

**Synonyms:** 4-bisphosphate 4-phosphatase; Inositol polyphosphate 4 phosphatase type II 105kDa; Inositol polyphosphate 4-phosphatase type II; Inpp4b; MGC132014; Type II inositol 3; 4 bisphosphate 4 phosphatase; Type II inositol-3

**UNIPROT ID:** O15327

**Background:** Catalyzes the hydrolysis of the 4-position phosphate of phosphatidylinositol 3,4-bisphosphate, inositol 1,3,4-trisphosphate and inositol 1,4-bisphosphate.

**Immunogen:** A synthesized peptide derived from human INPP4B

**Applications:** WB,IHC-P,FC

**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 FC: 1/50-1/100

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Clone ID:** -

**MW:** Calculated MW: 105 kDa; Observed MW: 105 kDa

**Isotype:** IgG

**Purification:** Affinity Chromatography

**Species Reactivity:** Human,Mouse,Rat

**Conjugation:** Unconjugated

**Modification:** Unmodified

**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

**Research Areas:** Signal Transduction

**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of INPP4B in MCF-7 lysates using INPP4B antibody.