

IQCD RABBIT PAB

Cat.#: S219073

Product Name: Anti-IQCD Rabbit Polyclonal Antibody

Synonyms: DRC10; CFAP84; 4933433C09Rik

UNIPROT ID: Q96DY2 (Gene Accession - BC013151)

Background: IQCD contains 1 IQ domain. The function of the IQCD protein is not known.

Immunogen: Fusion protein of human IQCD

Applications: ELISA, IHC

Recommended Dilutions: IHC: 50-300; ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human colorectal cancer tissue using 219073 (IQCD Antibody) at a dilution of 1/90 (Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human colorectal cancer tissue is first treated with the fusion protein and then with 219073 (Anti-IQCD Antibody) at dilution 1/90.