




KAT1 RABBIT MAB**Cat.#:** N262440**Product Name:** Anti-KAT1 Rabbit Monoclonal Antibody**Synonyms:** hat1; KAT1**UNIPROT ID:** O14929**Background:** Acetylates soluble but not nucleosomal histone H4 at 'Lys-5' (H4K5ac) and 'Lys-12' (H4K12ac) and, to a lesser extent, acetylates histone H2A at 'Lys-5' (H2AK5ac). Has intrinsic substrate specificity that modifies lysine in recognition sequence GXGKXG.**Immunogen:** Recombinant protein of human KAT1**Applications:** WB,IHC-F,IHC-P,ICC/IF,IP**Recommended Dilutions:** WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 IP: 1/20**Host Species:** Rabbit**Clonality:** Rabbit Monoclonal**Clone ID:** R02-4J1**MW:** Calculated MW: 50 kDa; Observed MW: 45 kDa**Isotype:** IgG**Purification:** Affinity Purified**Species Reactivity:** Human,Mouse,Rat**Conjugation:** Unconjugated**Modification:** Unmodified**Constituents:** PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide**Research Areas:** Epigenetics and Nuclear Signaling**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing



Immunocytochemistry analysis of KAT1 (green) in Hela using KAT1 antibody, and DAPI (blue)



Western blot analysis of KAT1 in K562, C6, 3T3, Hela lysates using KAT1 antibody.



Immunohistochemistry analysis of paraffin-embedded Human tonsil using KAT1/HAT1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.