

KCTD7 RABBIT PAB

Cat.#: N225227

Product Name: Anti-KCTD7 Rabbit pAb

Synonyms: BTB/POZ domain-containing protein KCTD7

UNIPROT ID: Q96MP8

Background: May be involved in the control of excitability of cortical neurons.

Immunogen: Synthesized peptide derived from human KCTD7. at AA range: 181-230

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 33 kDa; Observed MW: 33 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of KCTD7 in mouse brain , mouse spleen lysates using KCTD7 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
