

LCAT RABBIT PAB

Cat.#: S216613

Product Name: Anti-LCAT Rabbit Polyclonal Antibody

Synonyms:

UNIPROT ID: P04180 (Gene Accession - BC014781)

Background: This gene encodes the extracellular cholesterol esterifying enzyme, lecithin-cholesterol acyltransferase. The esterification of cholesterol is required for cholesterol transport. Mutations in this gene have been found to cause fish-eye disease as well as LCAT deficiency.

Immunogen: Fusion protein of human LCAT

Applications: ELISA, IHC

Recommended Dilutions: IHC: 25-100; ELISA: 2000-5000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Metabolism, Cancer, Cardiovascular

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human liver cancer tissue using 216613(LCAT Antibody) at a dilution of 1/30 (Secreted).



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
