

LIV-1(DMC488) IGG1 CHIMERIC MAB

目录: 28561

产品名称: LIV-1(DMC488) IgG1 Chimeric Monoclonal Antibody

基因符号: LIV1, SLC39A6, ZIP-6

描述: LIV-1 antibody(DMC488) IgG1 Chimeric Monoclonal Antibody

背景: Zinc is an essential cofactor for hundreds of enzymes. It is involved in protein, nucleic acid, carbohydrate, and lipid metabolism, as well as in the control of gene transcription, growth, development, and differentiation. SLC39A6 belongs to a subfamily of proteins that show structural characteristics of zinc transporters (Taylor and Nicholson, 2003 [PubMed 12659941]).[supplied by OMIM, Mar 2008]

经过测试的应用: Flow Cyt

推荐稀释度: Flow Cyt 1/100

种属反应性: Rabbit

亚型: Rabbit/Human Fc chimeric IgG1

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: Human LIV-1

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).




Figure 1. Flow cytometry analysis with 1 µg/mL LIV-1 (DMC488) mAb on RPMI 8226 cell line.

