

LLGL1 RABBIT MAB

Cat.#: N263407

Product Name: Anti-LLGL1 Rabbit Monoclonal Antibody

Synonyms: DLG4; HUGL; LLGL; Lgll; MglI; HUGL1; HUGL-1

UNIPROT ID: Q15334

Background: Cortical cytoskeleton protein found in a complex involved in maintaining cell polarity and epithelial integrity.

Immunogen: Recombinant protein of human LLGL1

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-0B8

MW: Calculated MW: 115 kDa; Observed MW: 115 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Neuroscience

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of LLGL1 in A549, U2OS lysates using LLGL1 antibody.