

MAP3K11 RABBIT PAB

Cat.#: S220716

Product Name: Anti-MAP3K11 Rabbit Polyclonal Antibody

Synonyms: MLK3; PTK1; SPRK; MLK-3; MEKK11

UNIPROT ID: Q16584 (Gene Accession - NP_002410)

Background: The protein encoded by this gene is a member of the serine/threonine kinase family. This kinase contains a SH3 domain and a leucine zipper-basic motif. This kinase preferentially activates MAPK8/JNK kinase, and functions as a positive regulator of JNK signaling pathway. This kinase can directly phosphorylate, and activates I κ B kinase alpha and beta, and is found to be involved in the transcription activity of NF- κ B mediated by Rho family GTPases and CDC42.

Immunogen: Synthetic peptide of human MAP3K11

Applications: ELISA, WB, IHC

Recommended Dilutions: IHC: 50-200;WB: 1000-5000;ELISA: 5000-10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Isotype: Immunogen-specific rabbit IgG


Purification: Antigen affinity purification

Species Reactivity: Human, Mouse, Rat


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

Research Areas: Signal Transduction


Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin embedded Human liver cancer tissue using 220716 (MAP3K11 Antibody) at a dilution of 1/50 (Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human liver cancer tissue is first treated with the synthetic peptide and then with 220716 (Anti-MAP3K11 Antibody) at dilution 1/50.



Gel: 6% SDS-PAGE, Lysate: 40 μ g;
Lane 1-5: K562, NIH/3T3, PC3, LO2, MCF7 cell lysates;
Primary antibody: 220716 (MAP3K11 Antibody) at dilution 1/800;
Secondary antibody: HRP-conjugated Goat anti rabbit IgG at 1/5000 dilution;
Exposure time: 20 seconds