

MCM3 RABBIT PAB

Cat.#: N225981

Product Name: Anti-MCM3 Rabbit pAb

Synonyms: MCM3; HCC5; P1-MCM3; P1.h; RLFB

UNIPROT ID: P25205

Background: Acts as component of the MCM2-7 complex (MCM complex) which is the putative replicative helicase essential for 'once per cell cycle' DNA replication initiation and elongation in eukaryotic cells.

Immunogen: A synthesized peptide derived from human MCM3

Applications: WB,IHC-P,ICC/IF,FC

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IF: 1/50-1/200 FC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 91 kDa; Observed MW: 100 kDa

Isotype: IgG

Purification: Affinity Chromatography

Species Reactivity: Human,Mouse,Rat

Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Epigenetics and Nuclear Signaling

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MCM3 in (1) K562 lysates; (2) RAW264.7 lysates using MCM3 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
