

MGST1 RABBIT MAB

Cat.#: N261580

Product Name: Anti-MGST1 Rabbit Monoclonal Antibody

Synonyms: Microsomal GST-I

UNIPROT ID: P10620

Background: Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Immunogen: A synthetic peptide of human Microsomal Glutathione S-transferase 1

Applications: WB,IHC-P

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R06-1E9

MW: Calculated MW: 18 kDa; Observed MW: 18 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Cell Biology

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human Cholangiocarcinoma using Microsomal Glutathione Stransferase 1 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of Microsomal GST1 in rat Liver lysates using MGST1 antibody.