

MGST3 RABBIT MAB

Cat.#: N262515

Product Name: Anti-MGST3 Rabbit Monoclonal Antibody

Synonyms: Glutathione peroxidase MGST3; LTC4 synthase MGST3; Microsomal GST-III

UNIPROT ID: O14880

Background: Catalyzes oxydation of hydroxy-fatty acids.

Immunogen: A synthetic peptide of human MGST3

Applications: WB,IHC-P,IP

Recommended Dilutions: WB: 1/500-1/1000 IHC: 1/50-1/100 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R07-6A3

MW: Calculated MW: 17 kDa; Observed MW: 17 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunohistochemistry analysis of paraffin-embedded Human colon cancer using MGST3 antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.



Western blot analysis of MGST3 in K562, HeLa lysates using MGST3 antibody.