

## MIA RABBIT PAB

**Cat.#:** S219202

**Product Name:** Anti-MIA Rabbit Polyclonal Antibody

**Synonyms:** CD-RAP

**UNIPROT ID:** Q16674 (Gene Accession - BC005910 )

**Background:** Elicits growth inhibition on melanoma cells in vitro as well as some other neuroectodermal tumors, including gliomas.

**Immunogen:** Fusion protein of human MIA

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 100-300; ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human, Mouse, Rat


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Signal Transduction, Cancer, Immunology


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing




Immunohistochemistry analysis of paraffin embedded Human cervical cancer tissue using 219202(MIA Antibody) at a dilution of 1/100(Secreted).



In comparison with the IHC on the left, the same paraffin-embedded Human cervical cancer tissue is first treated with the fusion protein and then with 219202(Anti-MIA Antibody) at dilution 1/100.



The image on the left is immunohistochemistry of paraffin-embedded Human colorectal cancer tissue using 219202(Anti-MIA Antibody) at a dilution of 1/100.



In comparison with the IHC on the left, the same paraffin-embedded Human colorectal cancer tissue is first treated with fusion protein and then with D226041(Anti-MIA Antibody) at dilution 1/100.