

## MID2 RABBIT PAB

**Cat.#:** S219518

**Product Name:** Anti-MID2 Rabbit Polyclonal Antibody

**Synonyms:** FXY2; RNF60; TRIM1; MRX101

**UNIPROT ID:** Q9UJV3 (Gene Accession - BC017707 )

**Background:** The protein encoded by this gene is a member of the tripartite motif (TRIM) family. The TRIM motif includes three zinc-binding domains, a RING, a B-box type 1 and a B-box type 2, and a coiled-coil region. The protein localizes to microtubular structures in the cytoplasm. Alternate splicing of this gene results in two transcript variants encoding different isoforms.

**Immunogen:** Fusion protein of human MID2

**Applications:** ELISA, IHC

**Recommended Dilutions:** IHC: 100-300; ELISA: 5000-10000

**Host Species:** Rabbit

**Clonality:** Rabbit Polyclonal

**Isotype:** Immunogen-specific rabbit IgG


**Purification:** Antigen affinity purification

**Species Reactivity:** Human, Mouse


**Constituents:** PBS (without Mg<sup>2+</sup> and Ca<sup>2+</sup>), pH 7.4, 150 mM NaCl, 0.05% Sodium Azide and 40% glycerol

**Research Areas:** Neuroscience, Signal Transduction, Epigenetics and Nuclear Signaling


**Storage & Shipping:** Store at -20°C. Avoid repeated freezing and thawing




Immunohistochemistry analysis of paraffin embedded Human thyroid cancer tissue using 219518 (MID2 Antibody) at a dilution of 1/95 (Cytoplasm).



In comparison with the IHC on the left, the same paraffin-embedded Human thyroid cancer tissue is first treated with the fusion protein and then with 219518 (Anti-MID2 Antibody) at dilution 1/95.



The image on the left is immunohistochemistry of paraffin-embedded Human gastric cancer tissue using 219518 (Anti-MID2 Antibody) at a dilution of 1/95.



In comparison with the IHC on the left, the same paraffin-embedded Human gastric cancer tissue is first treated with fusion protein and then with D227456 (Anti-MID2 Antibody) at dilution 1/95.