

MKNK2 RABBIT PAB

Cat.#: N225456

Product Name: Anti-MKNK2 Rabbit pAb

Synonyms: MKNK2; GPRK7; MNK2; MAP kinase-interacting serine/threonine-protein kinase 2; MAP kinase signal-integrating kinase 2; MAPK signal-integrating kinase 2; Mnk2

UNIPROT ID: Q9HBH9

Background: Serine/threonine-protein kinase that phosphorylates SFPQ/PSF, HNRNPA1 and EIF4E.

Immunogen: Synthesized peptide derived from the Internal region of human Mnk2.

Applications: WB,ELISA

Recommended Dilutions: WB: 1/500-1/1000 ELISA: 1/10000

Host Species: Rabbit

Clonality: Rabbit Polyclonal

Clone ID: -

MW: Calculated MW: 52 kDa; Observed MW: 52 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse,Rat


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MKNK2 in HeLa lysates using MKNK2 antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
