

## MMAE(8C4) RABBIT MAB

目录: 28571

产品名称: MMAE(8C4) Rabbit Monoclonal Antibody

基因符号: N.A.

描述: MMAE antibody(8C4) Rabbit Monoclonal Antibody

背景: N.A.

经过测试的应用: ELISA

推荐稀释度: ELISA 1/5000-10000

种属反应性: Rabbit

亚型: Rabbit IgG

纯化: Purified from cell culture supernatant by affinity chromatography

种属反应性: N.A. MMAE

成分: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

储存和运输: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

ELISA assay to evaluate Anti-MMAE antibody  
0.2µg Human IgG-MMAE per well




Figure1. Elisa plates were pre-coated with IgG-MMAE (0.2µg/per well). Serial diluted MMAE monoclonal antibody (28571) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC50 for MMAE monoclonal antibody binding with IgG-MMAE is 1.356ng/ml.

