

MMAE(8C4) RABBIT MAB

Cat.#: 28571

Product Name: Anti-MMAE(8C4) Rabbit Monoclonal Antibody

Synonyms: N.A.

Description: Anti-MMAE antibody(8C4) Rabbit Monoclonal Antibody

Background: N.A.

Applications: ELISA

Recommended Dilutions: ELISA 1/5000-10000

Host Species: Rabbit

Isotype: Rabbit IgG

Purification: Purified from cell culture supernatant by affinity chromatography

Species Reactivity: N.A. MMAE

Constituents: Lyophilized from sterile PBS, pH 7.4. 5 % – 8% trehalose is added as protectants before lyophilization.

Storage & Shipping: Store at -20°C to -80°C for 12 months in lyophilized form. After reconstitution, if not intended for use within a month, aliquot and store at -80°C (Avoid repeated freezing and thawing).

ELISA assay to evaluate Anti-MMAE antibody
0.2 μg Human IgG-MMAE per well




Figure1. Elisa plates were pre-coated with IgG-MMAE (0.2 μg /per well). Serial diluted anti-MMAE monoclonal antibody (28571) solutions were added, washed, and incubated with secondary antibody before Elisa reading. From above data, the EC₅₀ for anti-MMAE monoclonal antibody binding with IgG-MMAE is 1.356ng/ml.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
