

MVK RABBIT MAB

Cat.#: N263442

Product Name: Anti-MVK Rabbit Monoclonal Antibody

Synonyms: MVK; Mevalonate kinase; MK

UNIPROT ID: Q03426

Background: Catalyzes the phosphorylation of mevalonate to mevalonate 5-phosphate, a key step in isoprenoid and cholesterol biosynthesis.

Immunogen: A synthetic peptide of human MVK

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R08-6D3

MW: Calculated MW: 42 kDa; Observed MW: 42 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: -

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of MVK in HeLa, A549, HL-60, U251, U87-MG lysates using MVK antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
