

MYC TAG MAB

MYC Tag 小鼠单克隆抗体

货号: 28001

规格: 100 μ L

Symbol: C-Myc

描述: MYC Tag 小鼠单克隆抗体

背景: The MYC antibody recognizes any protein containing Myc Tag.

免疫原: CEQKLISEEDL

经过测试的应用: ELISA, WB, IF, IHC

推荐稀释度:

ELISA 1:2000-1:5000

WB 1:1000-1:5000

IF 1:100-1:200

浓度: 1 mg/ml

种属反应性: Mouse

形式: Liquid

克隆性: Monoclonal

亚型: IgG

纯化:: Purified from ascites


Preservative: No

成分: PBS (without Mg^{2+} and Ca^{2+}), pH 7.4, 150 mM NaCl, 50% glycerol

种属反应性: Recognizes protein contain Myc .

储存条件:: Store at $-20^{\circ}C$. Avoid repeated freezing and thawing

Western blot:



Western blot: Recombinant c-Myc tag fusion protein was diluted with c-Myc mouse monoclonal antibody (1:1000), and the protein loading amount was 50 ng, and Western blot analysis was performed.

