

MYH1 (9C7) MOUSE MAB

Cat.#: N261219

Product Name: Anti-MYH1 (9C7) Mouse Monoclonal Antibody

Synonyms: fast skeletal muscle; MYH3; HEMHC; Myosin-1; Myosin heavy chain 2x; MyHC-2x; Myosin heavy chain IIx/d; MyHC-IIx/d; Myosin heavy chain; skeletal muscle; adult 1

UNIPROT ID: P12882

Background: Muscle contraction.

Immunogen: Synthetic Peptide of Myosin Heavy Chain

Applications: ICC/IF

Recommended Dilutions: IF: 1/50-1/200

Host Species: Mouse

Clonality: Mouse Monoclonal

Clone ID: 9C7-4A6-5H10

MW: -

Isotype: IgG1

Purification: Affinity Purified

Species Reactivity: Human, Mouse, Rat, Fruitfly, Nematode


Conjugation: Unconjugated

Modification: Unmodified

Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Immunofluorescence analysis of MYH1 (9C7) in Immunofluorescence Staining of nematode tissue using MYH1 (9C7) antibody.



Product Description

Pioneering GTPase and Oncogene Product Development since 2010
