

MYH10 RABBIT MAB

Cat.#: N263458

Product Name: Anti-MYH10 Rabbit Monoclonal Antibody

Synonyms: NMMHCB; NMMHC-IIB

UNIPROT ID: P35580

Background: Cellular myosin that appears to play a role in cytokinesis, cell shape, and specialized functions such as secretion and capping.

Immunogen: A synthetic peptide of human non-muscle Myosin IIB/MYH10

Applications: WB,IP

Recommended Dilutions: WB: 1/500-1/1000 IP: 1/20

Host Species: Rabbit

Clonality: Rabbit Monoclonal

Clone ID: R05-6D8

MW: Calculated MW: 229 kDa; Observed MW: 229 kDa

Isotype: IgG

Purification: Affinity Purified

Species Reactivity: Human,Mouse


Conjugation: Unconjugated

Modification: Unmodified


Constituents: PBS (without Mg²⁺ and Ca²⁺), pH 7.3 containing 50% glycerol, 0.5% BSA and 0.02% sodium azide

Research Areas: Signal Transduction

Storage & Shipping: Store at -20°C. Avoid repeated freezing and thawing



Western blot analysis of non-muscle Myosin IIB/MYH10 in A549, HL-60, U251, U87-MG lysates using non-muscle Myosin IIB/MYH10 antibody.



Western blot analysis of non-muscle Myosin IIB/MYH10 in 3T3 lysates using MYH10 antibody.